

Rhodamine 6G*

Synonym: 2-[6-(ethylamino)-3-(ethylimino)-2,7-dimethyl-3H-xanthen-9-yl]-benzoic acid, ethyl ester, chloride or perchlorate; Basic Red 1 (Rhodamine 6G = ethyl ester; Rhodamine 590 = methyl ester; to the best of our knowledge, the two esters have equivalent performance)

Catalog No.: 05906 (chloride); 05905 (perchlorate)

CAS No.: 989-38-8 (05906); 13161-28-9 (05905)


MW: 479.02 (05906); 543.01 (05905)

Chemical Formula: C₂₈H₃₁ClN₂O₃ (05906); C₂₈H₃₁N₂O₃.ClO₄ (05905)

Appearance: Red solid (05906); red to violet (05905)

Molar Absorptivity (in ethanol): 9.70 x 10⁴ L mole⁻¹ cm⁻¹ (chloride form)

Structure:



Laser Dye Catalog No.	Lasing Wavelength		Pump Source (nm)	Solvent	Concentration (molar)	Abs λ-max	FI λ-max
	Max. (nm)	Range (nm)					
Rhodamine 6G (Available as the Chloride, Cl, 05906 ; the Perchlorate, ClO ₄ , 05905)	578	565-612	FL ³	Methanol	5 x 10 ⁻⁵	530 ^e	556 ^e
	584	570-618	FL ³	Ethanol	5 x 10 ⁻⁵		
	585	562-622	FL ⁶⁹	Methanol	4 x 10 ⁻⁵		
	586	563-625	FL ¹¹	Methanol	5 x 10 ⁻⁵		
	590		FL ⁶³	Methanol	8 x 10 ⁻⁵		
	596	577-614	FL ⁶⁹	MeOH/H ₂ O,1/3			
	598	577-625	FL ¹²	MeOH/H ₂ O,1/1	1.3 x 10 ⁻⁴		
		590-610	FL ¹⁸⁸	β-cyclodextrin/H ₂ O	2.3 x 10 ⁻⁴		
	600		FL ⁶³	4% LO/H ₂ O	1.2 x 10 ⁻⁴		
	610	585-633	FL ¹²	4% LO/H ₂ O	1.3 x 10 ⁻⁴		
	606		FL (Triaxial) ²²⁷	Acrylic Copolymer	1 x 10 ⁻⁴		
	580		KrF(248) ⁴⁴	Ethanol	1 x 10 ⁻³		
	590		KrF(248) ⁴⁶	p-Dioxane			
	574	563-615	XeCl(308) ¹¹⁴	Methanol	1.5 x 10 ⁻³		
	580	567-610	XeCl(308) ²⁰⁴	Ethanol	2.5 x 10 ⁻³ (osc), 1.8 x 10 ⁻³ (amp)		
	582	570-616	XeCl(308) ¹¹⁸	Ethanol	2.5 x 10 ⁻³ (osc), 3.8 x 10 ⁻⁵ (KR620)(amp)		
	583	566-610	XeCl(308) ¹¹⁰	Methanol	1.5 x 10 ⁻³		
	585	570-602	XeCl(308) ¹¹⁰	Methanol	1 x 10 ⁻³		
	591		XeCl(308) ¹¹²	Ethanol	4 x 10 ⁻³		
	586	570-614	XeF(351) ¹⁵⁴	Ethanol	5 x 10 ⁻³		
550		Nd:YAG(532) ⁵⁴	Methanol	3 x 10 ⁻⁴			
560	552-580	Nd:YAG(532) ⁵⁷	Methanol	2.2 x 10 ⁻⁴ (osc), 3.2 x 10 ⁻⁵ (amp)			
562	546-592	Nd:YAG(532) ⁵⁵	Methanol				


Laser Dye Catalog No.	Lasing Wavelength		Pump Source (nm)	Solvent	Concentration (molar)	Abs λ -max	Fl λ -max
	Max. (nm)	Range (nm)					
563	550-590		Nd:YAG(532) ⁵⁸	Methanol			
563	552-584		Nd:YAG(532) ⁵³	Methanol	120.6mg/l(osc), 51mg/l(amp)		
564			Nd:YAG(532) ⁵	Ethanol	3.7 x 10 ⁻⁴ (osc), 3 x 10 ⁻⁵ (amp)		
566	556-580		Nd:YAG(532) ¹¹⁰	Methanol	1 x 10 ⁻⁴		
566	559-576		Nd:YAG(532) ²³⁹	Ethanol	1.9 x 10 ⁻⁴		
574	563-597		Nd:YAG(355) ¹¹⁰	Methanol	6 x 10 ⁻⁴		
574	563-597		Nd:YAG(355) ²³⁹	Ethanol	8.4 x 10 ⁻⁴		
575	556-620		Nd:YAG(d,m-l,QS) ¹⁶⁸	EG			
575	565-600		Nd:YAG(532) ¹¹⁶	Ethanol	5 x 10 ⁻⁴		
577	567-602		Nd:YAG(355) ¹⁰⁹	Ethanol	2.5 x 10 ⁻³		
577			Nd:YAG(532,25kHz, 40 watts) ²²⁹	MeOH/H ₂ O, 1/1	1.75 x 10 ⁻⁴		
578	566-600		Nd:YAG(532) ¹¹⁰	Methanol	1 x 10 ⁻⁴ (R590) 1.3 x 10 ⁻⁵ (R610)		
567	557-590		N ₂ (337) ¹²²		+C540A		
576	555-618		N ₂ (337) ⁷³	Ethanol			
579	568-605		N ₂ (337) ⁵	Ethanol	5 x 10 ⁻³		
585	571-616		N ₂ (337) ¹¹⁴	Ethanol	4.2 x 10 ⁻³		
596	569-635		N ₂ (337) ⁹⁰	Ethanol	5.3 x 10 ⁻³		
596	575-625		N ₂ (337) ¹⁸³	Methanol	55mg/20ml		
590	570-650		Ar(458,514) ¹⁷	EG	2 x 10 ⁻³		
593	573-640		Ar(Blue/Green) ¹²³	LO/EG,3/1 + COT			
598	566-640		Ar(458-514) ²⁰⁶	EG/MeOH,9.5/0.5	2.1 x 10 ^{-3**}		
600	567-657		Ar(cw) ¹⁴	EG			
602	560-654		Kr(Blue/Green) ⁶⁸	MeOH/EG	80% pump absorption		
567	555-584		Cu(511) ¹⁵³	Methanol	4 x 10 ⁻⁴		
572	599-606		Cu(511) ¹⁵³	Methanol	9.6 x 10 ⁻⁴		
572	-564-600-		Cu(511,578) ²⁸	Ethanol	1 x 10 ⁻³		
585	563-607		Cu(511) ¹⁷⁵	Methanol	4 x 10 ⁻⁴		
590	575-614		Cu(511) ¹⁵³	Methanol	8.8 x 10 ⁻⁴ (R590)+ 2.1 x 10 ⁻⁴ (KR620)		

** This represents a maximum value. Concentration should be adjusted to 80-85% absorption of the pump light.


E = Ethanol; EG = Ethylene Glycol; MeOH = Methanol; LO = Ammonyx LO; COT = Cyclooctatetraene

* Equivalent species may be provided or substituted

Rhodamine 6G Chloride in Methanol



Rhodamine 6G Perchlorate in Methanol



The information presented above is believed to be accurate but is not a specification. The customer is fully responsible for determining the suitability of this product for use in their application. Exciton, Inc. does not represent that the information is sufficient or complete for any specific application.

REFERENCES:

3. Phase-R Corporation, Box G-2 Old Bay Rd., New Durham, NH 03855
5. Laser Photonics, Inc., 12351 Research Parkway, Orlando, FL 32826, formerly, Molelectron Corporation and Cooper LaserSonics, Inc.
11. Lasing Characteristics of Seventeen Visible-Wavelength Dyes using a Coaxial-Flashlamp-Pumped Laser, J.B. Marling, J.H. Hawley, E.M. Liston and W.B. Grant, *Appl. Optics*, 13(10), 2317 (1974). a. With Rhodamine 6G
12. Chromatix, 560 Oak Meade Parkway, Sunnyvale, CA 94086
14. CW Laser Emission Spanning the Visible Spectrum, J.M. Yarborough, *Appl. Phys. Lett.*, 24(12), 629 (1974). a. With Rhodamine 6G
17. Spectra-Physics, 1250 W. Middlefield Road, Mountain View, CA 94039
28. Efficient, High Average Power Dye Amplifiers Pumped by Copper Vapor Lasers, R.S. Hargrove and T. Kan, *IEEE J. Quantum Electron.*, QE13, 28D (1977)
44. Some Characteristics of Efficient Dye Laser Emission Obtained By Pumping at 248 nm with a High-Power KrF* Discharge Laser, V.I. Tomin, A.J. Alcock, W.J. Sarjeant, and K.E. Leopold, *Optics Commun.*, 26(3), 396 (1978)
46. Characterization of Dye Laser Pumping Using a High-Power KrF Excimer Laser at 248 nm, T.J. McKee, and D.J. James, to be published September 1979 in *Canadian J. Physics*
53. Continuum, 3150 Central Expressway, Santa Clara, CA 95051, formerly, Quantel International
54. W. R. Green, private commun., 1977
55. A High-Power Dye-Laser Pumped by the Second Harmonic of a Nd-YAG Laser, W. Hartig, *Optics Commun.*, 27(3), 447 (1978)
57. Quanta-Ray, Note: Quanta-Ray is now incorporated as a part of Spectra-Physics, 1250 W. Middlefield Road, Mountain View, CA 94039
58. J.K. Lasers Ltd., Somers Road, Rugby, Warwickshire, U. K.
63. High Energy Pulsed Dye Lasers for Atmospheric Sounding, J.Y. Allain, *Appl. Optics*, 18(3), 287 (1979)
68. Coherent Inc., 3210 Porter Dr., Palo Alto, CA 94304
69. Candela Laser Corporation, 530 Boston Post Road, Wayland, MA 01778-1833
73. Laser Dye DCM, Spectral Properties, Synthesis and Comparison with other Dyes in the Red, P.R. Hammond, *Optics Commun.*, 29(3), 331 (1979)
90. Jobin Yvon, 16-18 rue du Canal B.P. 118, 91163 Longjumeau Cedex France
109. Tuning Ranges of 355 nm Pumped Dyes from 410-715 nm, D.M. Guthals and J.W. Nibler, *Optics Commun.*, 29(3), 322 (1979)
110. Lumonics Inc., 105 Schneider Road, Kanata, (Ottawa), Ontario, Canada K2K 1Y3
112. Efficient Dye Lasers Pumped by an XeCl Excimer Laser, O. Uchino, T. Mizunami, M. Maeda and Y. Miyazoe, *Appl. Phys.*, 19, 35 (1979)
114. Optimization of Spectral Coverage in an Eight-Cell Oscillator-Amplifier Dye Laser Pumped at 308nm, F. Bos, *Appl. Optics*, 20, 3553 (1981)
116. Versatile High-Power Single-Longitudinal-Mode Pulsed Dye Laser, F. Bos, *Appl. Optics*, 20(10), 1886 (1981)
118. The XeCl Excimer Laser: A Powerful and Efficient UV Pumping Source for Tunable Dye Lasers, H. Telle, W. Huffer and D. Basting, *Optics Commun.*, 38(5,6), 402 (1981)
122. Increased Gain Through Identification and Alleviation of Dye Self-Absorption in Laser-Pumped Dye Laser, R.E. Drullinger, *Optics Commun.*, 39(4), 263 (1981)
123. Powerful Single-Frequency Ring Dye Laser Spanning the Visible Spectrum, T.F. Johnston, Jr., R.H. Brady and W. Proffitt, *Appl. Optics*, 21(13), 2307 (1982)
153. Cooper LaserSonics, Inc. 5674 Sonoma Drive, Pleasanton, CA 94566
154. Dye Laser Radiation in the 370-760nm Region Pumped by a XeF Excimer Laser, T.C. Eschrich and T.J. Morgan, *Applied Optics*, 24(7), 937 (1985)
168. Simple Picosecond Dye Laser System Pumped by a Frequency Doubled, Optically Compressed Q-Switched Mode-Locked Nd:YAG Laser, A.S.L. Gomes and J.R. Taylor, *Optics Commun.*, 55(6), 435 (1985)
175. CVL-Pumped Dye Laser For Spectroscopic Application, M. Broyer, J. Chevaleyre, G. Delacretaz and L. Wöste, *App. Phys. B*, 35, 31 (1984)
183. Laser Science, Inc., 26 Landsdowne Street, Cambridge, MA 02139
188. Effect of Cyclodextrins and Water-Soluble Polymers on the Spectroscopic and Lasing Properties of Rhodamine 6G, M.M. Asimov, V.P. Chuev, S.N. Kovalenko, V.M. Nikitchenko and A.N. Rubinov, *Opt. Spectrosc. (USSR)*, 70(3), 320 (1991)
204. Questek, Inc., 44 Manning Road, Billerica, MA 01821 (Tuning Curves for Model 5200B Dye Laser, PDL-3)
206. Coherent Inc., 3210 Porter Dr., Palo Alto, CA 94304; (599 Composite Tuning Curves, 1992; The concentration shown represents a maximum value. The final concentration should be adjusted to obtain 80-85% absorption of the pump light.)
227. Progress in Solid State Dye Laser Development, R.E. Hermes, Proceedings of the Int. Conf. on Lasers '90, STS Press, (1991)
229. Continuous-wave Dye Lasers in the DCM Region, P. Hammond and D. Cooke, *Appl. Optics* 31(33), 7095 (1992)



2150 Bixby Road
Lockbourne, OH 43137
Tel: 614.492.5610
E-mail: info.exciton@luxotticaretail.com
www.exciton.luxottica.com

239. P. Jauernik, private commun., Sirah Laser- und Plasmatechnik, 2003.

For a current list of biology, biological stain, or biochemistry references for Rhodamine 590 from PubMed, click on the following link:
[Rhodamine 590 or Rhodamine 6G](#)