

OXAZINE 725 PERCHLORATE or CHLORIDE

Synonym: 3,7-bis(diethylamino)phenoxazin-5-ium perchlorate; Oxazine 1
Phenoxazin-5-ium, 3,7-bis(diethylamino)-, perchlorate (1:1)
Phenoxazin-5-ium, 3,7-bis(diethylamino)-, chloride (1:1)

Catalog No.: 07250 (perchlorate), 07251 (chloride)

CAS No.: 24796-94-9 (perchlorate); 33203-82-6 (chloride); 47367-75-9 (cation only)


Chemical Formula: C₂₀H₂₆N₃O.ClO₄

C₂₀H₂₆N₃O.Cl

MW: 423.90 (07250), 359.90 (07251)

Appearance: Green bronzy crystals

Structure:




Lasering Wavelength


Max. (nm)	Range (nm)	Pump Source (nm)	Solvent	Concentration (molar)	Abs λ-max	FI λ-max
681		FL→R610(622) ⁶⁵	CH ₂ Cl ₂	4 x 10 ⁻⁵	645 ^e	680 ^e
695		FL→R610(622) ⁶⁵	DMSO	1 x 10 ⁻⁴		
715		FL ²⁷	Ethanol			
725	705-745	FL ¹¹	MeOH/R590			
740	720-758	FL ³	CH ₂ Cl ₂	3.3 x 10 ⁻⁵		
725	690-765	XeCl(308) ¹¹⁴	Methanol	2.4 x 10 ⁻³		
730	692-751	XeCl(308) ¹¹⁰	Methanol	1 x 10 ⁻³		
734	692-768	XeCl(308) ¹¹⁸	Ethanol	2.1 x 10 ⁻³		
742	702-772	XeCl(308) ¹¹⁰	Methanol	2 x 10 ⁻³		
735	705-758	XeF(351) ¹⁵⁴	Ethanol	5 x 10 ⁻³ (OX725)+ 5 x 10 ⁻³ (R610)		
675	657-695	Nd:YAG(532) ¹¹⁶	Methanol	4.9 x 10 ⁻⁵ (OX725), 9.5 x 10 ⁻⁴ (R640)		
690		Nd:YAG(532) ³³	CH ₂ Cl ₂	4 x 10 ⁻⁴		
690	671-712	Nd:YAG(532) ¹¹⁶	Methanol	2 x 10 ⁻⁵ (OX725), 9.5 x 10 ⁻⁴ (R640)		
720	695-755	Nd:YAG(532) ¹¹⁶	Methanol	1.2 x 10 ⁻³ (OX725), 6 x 10 ⁻⁴ (R640)		
724	695-761	Nd:YAG(532) ⁵³	Methanol			
725	705-750	N ₂ (337) ⁵	Ethanol	5 x 10 ⁻³ (R610), 5 x 10 ⁻³ (OX725)		
730	692-751	N ₂ (337) ¹¹¹	Ethanol/DMSO,96/4	2.3 x 10 ⁻³ (OX725), 2.5 x 10 ⁻³ (R610)		
730	692-751	N ₂ (337) ¹⁸³	DMSO	20mg/20ml		
750	736-765	N ₂ (337) ¹¹¹	DMSO	2 x 10 ⁻²		
723	688-800	Kr(Red) ⁶⁸ EG				
	687-826	Kr(647) ⁷¹	DMSO/EG,1/3 + COT	1 x 10 ⁻³		
745	645-810	Kr(647,676) ^{36b}	DMSO/EG or G	1.1 x 10 ⁻³		
750	695-801	Kr(647,676) ¹⁷	EG/DMSO,84/16	6 x 10 ⁻⁴		
729		HeNe(633,20mW) ¹⁷⁰	EtOH(-35°C)	1.6 x 10 ⁻³		
760	744-776	AlGInP (laser diode, 674) ²⁰⁸	EG	1.81x10 ⁻³		
760	750-770	AlGInP (laser diode, 674) ^{208a}	EG	1.81x10 ⁻³		

DMSO = Dimethylsulfoxide; EG = Ethylene Glycol; EtOH = Ethanol; G = Glycerol; COT = Cyclooctatetraene;
 MeOH = Methanol; CH₂Cl₂ = Methylene Chloride; e = Ethanol

Oxazine 725 Perchlorate in Methanol



Oxazine 725 Chloride in Methanol



The information presented above is believed to be accurate but is not a specification. The customer is fully responsible for determining the suitability of this product for use in their application. Exciton, Inc. does not represent that the information is sufficient or complete for any specific application.

REFERENCES:

3. Phase-R Corporation, Box G-2 Old Bay Rd., New Durham, NH 03855
5. Laser Photonics, Inc., 12351 Research Parkway, Orlando, FL 32826, formerly, Molelectron Corporation and Cooper LaserSonics, Inc.
11. Lasing Characteristics of Seventeen Visible-Wavelength Dyes using a Coaxial-Flashlamp-Pumped Laser, J.B. Marling, J.H. Hawley, E.M. Liston and W.B. Grant, *Appl. Optics*, 13(10), 2317 (1974). **a.** With Rhodamine 6G
17. Spectra-Physics, 1250 W. Middlefield Road, Mountain View, CA 94039
27. What's Ahead in Laser Dyes, K.H. Drexhage, *Laser Focus*, 9(3), 35 (1973)
33. A High-Power Dye Laser at 6700-7700 Å, K. Kato, *Optics Commun.*, 19(1), 18 (1976)
36. **a.** Spectra-Physics Laser Review, 4(1), April 1977; **b.** High Power CW Dye Laser Emission in the Near IR from 685 nm to 965 nm, K.M. Romanek, O. Hildebrand and E.Gobel, *Optics Commun.*, 21(1), 16 (1977)
53. Continuum, 3150 Central Expressway, Santa Clara, CA 95051, formerly, Quantel International
65. High-Power TEM₀₀ Tunable Laser System, R. Mahon, T.J. McIlrath and D.W. Koopman, *Appl. Optics*, 18(6), 891 (1979)
68. Coherent Inc., 3210 Porter Dr., Palo Alto, CA 94304
71. Generation of Near-Infrared Picosecond Pulses by Mode Locked Synchronous Pumping of a Jet-Stream Dye Laser, J. Kuhl, R. Lambrich and D. Von der Linde, *Appl. Phys. Lett.*, 31(10), 657 (1977)
110. Lumonics Inc., 105 Schneider Road, Kanata, (Ottawa), Ontario, Canada K2K 1Y3
111. Lasing Properties of Several Near-IR Dyes for a Nitrogen Laser-Pumped Dye Laser with an Optical Amplifier, B.M. Pierce and R.R. Birge, *IEEE J. Quantum Electron.*, QE18, 1164 (1982)
114. Optimization of Spectral Coverage in an Eight-Cell Oscillator-Amplifier Dye Laser Pumped at 308nm, F. Bos, *Appl. Optics*, 20, 3553 (1981)
116. Versatile High-Power Single-Longitudinal-Mode Pulsed Dye Laser, F. Bos, *Appl. Optics*, 20(10), 1886 (1981)
118. The XeCl Excimer Laser: A Powerful and Efficient UV Pumping Source for Tunable Dye Lasers, H. Telle, W. Huffer and D. Basting, *Optics Commun.*, 38(5,6), 402 (1981)
154. Dye Laser Radiation in the 370-760nm Region Pumped by a XeF Excimer Laser, T.C. Eschrich and T.J. Morgan, *Applied Optics*, 24(7), 937 (1985)
170. Continuous Wave Dye Laser Pumped by a HeNe Laser, E. Thiel, C. Zander and K. H. Drexhage, *Optics Commun.*, 60(6), 396 (1986)
183. Thermo Laser Science, 26 Landsdowne Street, Cambridge, MA 02139
208. Near-IR Dye Laser for Diode-Pumped Operation, R. Scheps, *IEEE J. Quantum Electron.* 31(1), 126 (1995)

For a current list of biology, biological stain, or biochemistry references for Oxazine 725 Perchlorate from PubMed, click on the following link:

[Oxazine 725 or Oxazine 1](#) (all references are listed under Oxazine 1 as of May 2006)