

### LDS 820/821

**Synonym:** LDS 820 – 2-[6-[4-(dimethylamino)phenyl]-1,3,5-hexatrienyl]-3-ethyl-benzothiazolium perchlorate;  
LDS 821 - (2-(6-(p-dimethylaminophenyl)-2,4-neopentylene-1,3,5-hexatrienyl)-3-ethylbenzothiazolium perchlorate);  
Styryl 9M

**Catalog No.:** 08200 (LDS 820); 08210 (LDS 821)


**CAS No.:** 76433-25-5 (08200); Not Available for 08210

**MW:** 460.98 (08200); 529.09 (08210)

**Chemical Formula:** C<sub>23</sub>H<sub>25</sub>N<sub>2</sub>S.ClO<sub>4</sub> (08200); C<sub>28</sub>H<sub>33</sub>N<sub>2</sub>S.ClO<sub>4</sub> (08210)

**Appearance:** Dark green crystals (08200); green crystals (08210)

**Structure:**



Lasing Wavelength Max. (nm)	Range (nm)	Pump Source (nm)	Solvent	Concentration (molar)	Abs λ-max	FI λ-max	% CE*
-----------------------------	------------	------------------	---------	-----------------------	-----------	----------	-------

All of the data for LDS 820/LDS 821 is interchangeable. LDS 821 is the preferred dye. LDS 820 is no longer available.

LDS 820 (08200):

841	810-860	FL <sup>141</sup>	PC/EG,9/1	2.5 x 10 <sup>-5</sup>			
810	784-844	Nd:YAG(532) <sup>134</sup>	Ethanol				
814	787-847	Nd:YAG(532) <sup>110</sup>	Methanol	1.5 x 10 <sup>-4</sup>			
818	775-865	Nd:YAG(532) <sup>57</sup>	Methanol	4.1 x 10 <sup>-4</sup> (osc), 3.9 x 10 <sup>-5</sup> (amp)			
825	798-866	Nd:YAG(532) <sup>134</sup>	DMSO				
840	800-865	Nd:YAG(532) <sup>5</sup>	PC,3%PC/EtOH	1 x 10 <sup>-3</sup> (LDS 820)(osc), 16.8mg/l(LDS 867)(amp)			
	775-840	Ar(m-l,514) <sup>136</sup>	PC/EG,15/85	1.9 x 10 <sup>-3</sup>			
	781-840	Ar(m-l,514) <sup>136</sup>	PC/EG,15/85	1.9 x 10 <sup>-3</sup>			
	792-880	Ar(m-l,514) <sup>136</sup>	PC/EG,15/85	1.9 x 10 <sup>-3</sup>			
822	784-900	Ar <sup>127a</sup>	PC/EG,15/85	1.2 x 10 <sup>-3</sup>			
	790-913	Ar(m-l,514) <sup>136</sup>	PC/EG,15/85	1.9 x 10 <sup>-3</sup>			
845	780-960	Ar(458-514) <sup>17</sup>	PC/EG,15/85	2 x 10 <sup>-3</sup>			

LDS 821 (08210):

834	817-842	FL <sup>69</sup>	Methanol	8.7 x 10 <sup>-5</sup>	574 <sup>m</sup>	750 <sup>m</sup>	--
848	824-867	FL <sup>69</sup>	DMSO	5.0 x 10 <sup>-5</sup>			--
818	785-850	XeCl(308) <sup>110</sup>	Methanol	6 x 10 <sup>-4</sup>			--
843	807-900	XeCl(308) <sup>204</sup>	DMSO	2.14 x 10 <sup>-3</sup> (osc), 1.28 x 10 <sup>-3</sup> (amp)			7
	805-840	Nd:YAG(532,m-l,1mj,32ps) <sup>172</sup>	Methanol	3.5 x 10 <sup>-3</sup> (cavity) 4 x 10 <sup>-4</sup> (amp)			--
812	780-844	Nd:YAG(532) <sup>53</sup>	Methanol	1.9 x 10 <sup>-4</sup> (osc), 2.5 x 10 <sup>-5</sup> (amp)			8.2
812	785-855	Nd:YAG(532) <sup>230</sup>	Methanol	189.5mg/l(osc), 72.8mg/l(amp)			--
815	791-839	Nd:YAG(532) <sup>239</sup>	Ethanol	2.5 x 10 <sup>-4</sup>			


818	785-851	Nd:YAG, (side-p,532) <sup>57</sup>	Methanol	125mg/l(osc), 32.5mg/l(amp)	10
818	785-851	Nd:YAG, (end-p,532) <sup>57</sup>	Methanol	125 mg/l(osc), 17mg/l(amp)	7
821		Nd:YAG(532,m-l) <sup>160</sup>	DMSO/PC/EG	$2.3 \times 10^{-3}$	11
826	780-874	Nd:YAG(532) <sup>53</sup>	Methanol	71.5mg/l(LDS 821)+ 28.6mg/l(LDS 867)(osc), 10mg/l(LDS 821)+ 9.3mg/l(LDS 867)(amp)	--
839	814-862	Nd:YAG(532) <sup>239</sup>	DMSO	$2.5 \times 10^{-4}$	
821	802-852	N <sub>2</sub> (337) <sup>137</sup>	PC	$3 \times 10^{-3}$ (osc), $3 \times 10^{-3}$ (amp)	--
815	782-930	Ar(514,m-l) <sup>136</sup>	PC/EG,15/85	$1.85 \times 10^{-3}$	29
840	790-940	Ar(458-514) <sup>206</sup>	PC/EG,3/7	$3.9 \times 10^{-3}$ *	14
843	780-960	Ar(459-514) <sup>17</sup>	PC/EG,15/85	$2 \times 10^{-3}$	11.3
880	793-923	Kr(647) <sup>128</sup>	PC/EG,1/4	1.4g/l	--
815	793-845	Cu(511,578) <sup>175</sup>	Methanol	$1.3 \times 10^{-3}$	14.4

\* This represents a maximum value. Concentration should be adjusted to 80-85% absorption of the pump light.  
 DMSO = dimethylsulfoxide, EG = ethylene glycol, PC = propylene carbonate, m = methanol

**NOTES:**

CE = Conversion efficiency reported by the manufacturer or literature sources. See reference (numbers indicated under pump source column)  
 -- = not reported or not available

**LDS 821 in Methanol**



The information presented above is believed to be accurate but is not a specification. The customer is fully responsible for determining the suitability of this product for use in their application. Exciton, Inc. does not represent that the information is sufficient or complete for any specific application.

**REFERENCES:**

5. Laser Photonics, Inc., 12351 Research Parkway, Orlando, FL 32826, formerly, Molelectron Corporation and Cooper LaserSonics, Inc.
17. Spectra-Physics, 1250 W. Middlefield Road, Mountain View, CA 94039
53. Continuum, 3150 Central Expressway, Santa Clara, CA 95051, formerly, Quantel International
57. Quanta-Ray, Note: Quanta-Ray is now incorporated as a part of Spectra-Physics, 1250 W. Middlefield Road, Mountain View, CA 94039
69. Candela Laser Corporation, 530 Boston Post Road, Wayland, MA 01778-1833
110. Lumonics Inc., 105 Schneider Road, Kanata, (Ottawa), Ontario, Canada K2K 1Y3
127. **a.** Cw Operation of Laser Dyes Styryl-9 and Styryl-11, J. Hoffnagle, L. Ph. Roesch, N. Schlumpf and A. Weis, *Optics Commun.*, 42, 267 (1982);  
**b.** K. Kato, see Reference 5 in 127 **a** ; **c.** K. Kato, unpublished results
128. F. Schellenberg, private commun., 1982; Dye Laser with HITC Optics and IR-140 Output Coupler Pumped by 3.8W Krypton-ion Laser
134. D. Heiman, private commun., 1983
136. High Efficiency Picosecond Pulse Generation in the 675-930nm Region from a Dye Laser Synchronously Pumped by an Argon-Ion Laser, P. Bado, C. Dupuy, K.R. Wilson, R. Boggy, J. Bowen and S. Westra, *Optic Commun.*, 46(3,4), 241 (1983)
137. Nitrogen Pumped LDS-821, P.Klein, private commun., 1983 (Molelectron UV 14 N<sub>2</sub> Laser/DL 14 Dye Laser - 3.5 mj/pulse-7% efficiency)
141. Flashlamp-Pumped Styryl 9 Dye Laser, Cheng-Huei Lin and B. Marshall, *Appl. Optics*, 23(14), 2228 (1984)
160. G. Olbright, private commun., 1988; Generation of Tunable Near-Infrared Amplified Femtosecond Laser Pulses and Time-Correlated White-Light Continuum, G.R. Olbright and G.R. Hadley, *J. Opt. Soc. Am. B*, 6(7), 1363 (1989)
172. Tunable Near-Infrared Picosecond Pulses from a Short-Cavity Dye Laser, K. Bohnert, T.F. Boggess, K. Mansour, D. Maxson and A.L. Smirl, *IEEE J. Quantum Electron.*, QE-22(12), 2195 (1986)
175. CVL-Pumped Dye Laser For Spectroscopic Application, M. Broyer, J. Chevalere, G. Delacretaz and L. Wöste, *App. Phys. B*, 35, 31 (1984)
204. Questek, Inc., 44 Manning Road, Billerica, MA 01821 (Tuning Curves for Model 5200B Dye Laser, PDL-3)
206. Coherent Inc., 3210 Porter Dr., Palo Alto, CA 94304; (599 Composite Tuning Curves, 1992; The concentration shown represents a maximum value. The final concentration should be adjusted to obtain 80-85% absorption of the pump light.)
230. Generation of 1.30- to 1.55 $\mu$ m Tunable Radiation from First Stokes Raman Shifting in Hydrogen, K.W. Aniolek, D.L. Miller, N.P. Cernansky, and K.G. Owens, *Appl. Spec.* 51(6), 820(1997)
239. P. Jauernik, private commun., Sirah Laser- und Plasmatechnik, 2003.

For a current list of biology, biological stain, or biochemistry references for LDS 821 from PubMed, click on the following link:  
[LDS 820/821 or Styryl 9M](#) (zero references in PubMed as of May 2006)