

## RHODAMINE 640 PERCHLORATE

**Synonym:** 9-(2-carboxyphenyl)-2,3,6,7,12,13,16,17-octahydro-1H,5H,11H,15H,-xantheno[2,3,4-ij:5,6,7-ij']-diquinolizin-4-ium perchlorate; Rhodamine 101

**Catalog No.:** 06400

**CAS No.:** 72102-91-1


**MW:** 591.05

**Chemical Formula:** C<sub>32</sub>H<sub>31</sub>N<sub>2</sub>O<sub>3</sub>.ClO<sub>4</sub>

**Appearance:** Dark green crystals with bronze sheen

**Molar absorptivity (at 567nm):** 10.50 x 10<sup>4</sup> L mole<sup>-1</sup> cm<sup>-1</sup>

**Structure:**




Max. Lasing Wavelength (nm)	Range (nm)	Pump Source (nm)	Solvent	Concentration (molar)	Abs λ-max	FI λ-max
630		FL <sup>27</sup>	Ethanol(basic)		575(a) <sup>m</sup>	594 <sup>m</sup>
635	610-670	FL <sup>69</sup>	Methanol	4 x 10 <sup>-5</sup>		
640		FL <sup>27</sup>	Ethanol(acidic)			
642	627-657	FL <sup>69</sup>	Methanol			
643	623-657	FL <sup>3</sup>	Ethanol	1 x 10 <sup>-4</sup>		
650		FL <sup>63</sup>	Methanol	1.2 x 10 <sup>-4</sup>		
652	620-687	FL <sup>12</sup>	MeOH/H <sub>2</sub> O,3/2	1.1 x 10 <sup>-4</sup>		
618	608-668	XeCl(308) <sup>114</sup>	Ethanol	1.3 x 10 <sup>-3</sup>		
623	613-672	XeCl(308) <sup>118</sup>	Ethanol	1.2 x 10 <sup>-3</sup> (osc)		
625	610-673	XeCl(308) <sup>204</sup>	Methanol	1.25 x 10 <sup>-3</sup> (osc), 7 x 10 <sup>-4</sup> (amp)		
602	589-623	Nd:YAG(532) <sup>55</sup>				
603	598-626	Nd:YAG(532) <sup>57</sup>	Methanol	2.4 x 10 <sup>-4</sup> (osc), 3.2 x 10 <sup>-5</sup> (amp)		
605	594-629	Nd:YAG(532) <sup>53</sup>	Methanol	R640(3.6x10 <sup>-5</sup> ) + R610(7.9x10 <sup>-5</sup> )(osc), R640(1.8x10 <sup>-5</sup> ) + R610(3.9x10 <sup>-5</sup> )(amp)		
611		Nd:YAG(532) <sup>54</sup>	Methanol	5 x 10 <sup>-4</sup>		
612	598-640	Nd:YAG(532) <sup>58</sup>				
613	605-630	Nd:YAG(532) <sup>53</sup>	Methanol	3.6 x 10 <sup>-4</sup> (osc), 1.9 x 10 <sup>-4</sup> (amp)		
619	607-640	Nd:YAG(532) <sup>110</sup>	Methanol	2.5 x 10 <sup>-4</sup>		
620	608-668	Nd:YAG(532) <sup>116</sup>	Ethanol	5 x 10 <sup>-4</sup>		
624	614-662	Nd:YAG(532) <sup>239</sup>	Ethanol	4.2 x 10 <sup>-4</sup>		
630	621-674	Nd:YAG(355) <sup>239</sup>	Ethanol	8.5 x 10 <sup>-4</sup>		
650	620-680	Nd:YAG(355) <sup>109</sup>	MeOH/H <sub>2</sub> O,3/2	3.5 x 10 <sup>-3</sup>		
640	620-680	N <sub>2</sub> (337) <sup>30</sup>	Ethanol	5 x 10 <sup>-3</sup>		
644	620-673	N <sub>2</sub> (337) <sup>50</sup>	Ethanol	5.7 x 10 <sup>-3</sup>		
652	620-678	N <sub>2</sub> (337) <sup>114</sup>	Ethanol	5.1 x 10 <sup>-3</sup>		
653	625-680	N <sub>2</sub> (337) <sup>183</sup>	Methanol	60mg/20ml		

659	626-700	N <sub>2</sub> (337) <sup>90</sup>	Ethanol	5 x 10 <sup>-3</sup>
671	634-704	N <sub>2</sub> (337) <sup>73</sup>	DMSO + HCl	
645	620-690	Ar(458-514) <sup>17</sup>	EG	1.5 x 10 <sup>-3</sup> (R640), 1.5 x 10 <sup>-3</sup> (R590)
648	608-710	Ar or Kr(568) <sup>68</sup>	MeOH/EG, 1/7.5	80% pump absorption
652	624-675	Ar(vis) <sup>87</sup>	EG	2 x 10 <sup>-3</sup> (R640), 1 x 10 <sup>-3</sup> (R590)
616	605-633	Cu(511,578) <sup>153</sup>	Methanol	1 x 10 <sup>-3</sup>
630	607-659	Cu(511,578) <sup>175</sup>	Methanol	

DMSO = Dimethylsulfoxide; EG = Ethylene glycol; HCl = Hydrochloric acid; MeOH/H<sub>2</sub>O = Methanol/water; m = methanol

### Rhodamine 640 Perchlorate in Methanol



The information presented above is believed to be accurate but is not a specification. The customer is fully responsible for determining the suitability of this product for use in their application. Exciton, Inc. does not represent that the information is sufficient or complete for any specific application.

### Quantum Yields and Lifetimes

Absorbance (nm)	Emission (nm)	Quantum Yield (max = 1.0)	Solvent	Lifetime (ns)	References, Notes
		1	Ethanol	5.3	R-2(16)
	588	0.96	Ethanol	4.6	R-2
			Basic Ethanol		R-8
	585-665	0.913	Fluorescent quantum yield relative to QSH (ff 0.55) at 25 degrees C; $1.34 \times 10^{-7} M$		
		1	Ethanol		S-4
			Methanol	4.9	R-2

### REFERENCES:

3. Phase-R Corporation, Box G-2 Old Bay Rd., New Durham, NH 03855
12. Chromatix, 560 Oak Meade Parkway, Sunnyvale, CA 94086
17. Spectra-Physics, 1250 W. Middlefield Road, Mountain View, CA 94039
27. What's Ahead in Laser Dyes, K.H. Drexhage, *Laser Focus*, 9(3), 35 (1973)
30. S. Woodruff and D. Ahlgren, private commun., 1977
50. G. Holtom, private commun., 1978
53. Continuum, 3150 Central Expressway, Santa Clara, CA 95051, formerly, Quantel International
54. W. R. Green, private commun., 1977
55. A High-Power Dye-Laser Pumped by the Second Harmonic of a Nd-YAG Laser, W. Hartig, *Optics Commun.*, 27(3), 447 (1978)
57. Quanta-Ray, Note: Quanta-Ray is now incorporated as a part of Spectra-Physics, 1250 W. Middlefield Road, Mountain View, CA 94039
58. J.K. Lasers Ltd., Somers Road, Rugby, Warwickshire, U. K.
63. High Energy Pulsed Dye Lasers for Atmospheric Sounding, J.Y. Allain, *Appl. Optics*, 18(3), 287 (1979)
68. Coherent Inc., 3210 Porter Dr., Palo Alto, CA 94304
69. Candela Laser Corporation, 530 Boston Post Road, Wayland, MA 01778-1833
73. Laser Dye DCM, Spectral Properties, Synthesis and Comparison with other Dyes in the Red, P.R. Hammond, *Optics Commun.*, 29(3), 331 (1979)
87. The Selective Excitation of Lithium Isotopes by Intracavity Nonlinear Absorption in a CW Dye Laser, M. Yamashita, M. Kasamatsu, H. Kashiwagi and K. Mashida, *Optics Commun.*, 26(3), 343 (1978)
90. Jobin Yvon, 16-18 rue du Canal B.P. 118, 91163 Longjumeau Cedex France
109. Tuning Ranges of 355 nm Pumped Dyes from 410-715 nm, D.M. Guthals and J.W. Nibler, *Optics Commun.*, 29(3), 322 (1979)
114. Optimization of Spectral Coverage in an Eight-Cell Oscillator-Amplifier Dye Laser Pumped at 308nm, F. Bos, *Appl. Optics*, 20, 3553 (1981)
116. Versatile High-Power Single-Longitudinal-Mode Pulsed Dye Laser, F. Bos, *Appl. Optics*, 20(10), 1886 (1981)
118. The XeCl Excimer Laser: A Powerful and Efficient UV Pumping Source for Tunable Dye Lasers, H. Telle, W. Huffer and D. Basting, *Optics Commun.*, 38(5,6), 402 (1981)

153. Cooper LaserSonics, Inc. 5674 Sonoma Drive, Pleasanton, CA 94566
175. CVL-Pumped Dye Laser For Spectroscopic Application, M. Broyer, J. Chevalayre, G. Delacretaz and L. Wöste, *App. Phys. B*, 35, 31 (1984)
183. Laser Science, Inc., 26 Landsdowne Street, Cambridge, MA 02139
204. Questek, Inc., 44 Manning Road, Billerica, MA 01821 (Tuning Curves for Model 5200B Dye Laser, PDL-3)
239. P. Jauernik, private commun., Sirah Laser- und Plasmatechnik, 2003.
- R-2. Photophysical Properties of Laser Dyes: Picosecond Laser Flash Photolysis Studies of Rhodamine 6G, Rhodamine B and Rhodamine 101, P.C. Beaumont D.G. Johnson, and B.J. Parsons, *J. of Chem. Soc., Farady Trans 89(23)*, 4185 (1993), <https://doi.org/10.1039/FT9938904185>
- R-8. Fluorescence Quantum Yields of Some Rhodamine Dyes, R.F. Kubin and A.N. Fletcher, *J. of Luminescence* 27, 455 (1982), [https://doi.org/10.1016/0022-2313\(82\)90045-X](https://doi.org/10.1016/0022-2313(82)90045-X)
- S-4. Fluorescence Quantum Yields of a Series of Red and Near-Infrared Dyes Emitting at 600-1000nm, K. Rurack and M. Spieles, *Anal. Chem.*, 83, 1232 (2011), <https://doi.org/10.1021/ac101329h>

For a current list of biology, biological stain, or biochemistry references for Rhodamine 640 from PubMed, click on the following link:

[Rhodamine 640 or Rhodamine 101](#) (all references are listed under Rhodamine 101 in PubMed as of May 2006)